

For Research Use Only

ANXA13 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 85955-3-PBS



Basic Information

Catalog Number: 85955-3-PBS	GenBank Accession Number: BC125158	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 312	CloneNo.: 250350F2
Isotype: IgG	UNIPROT ID: P27216	
Immunogen Catalog Number: AG18581	Full Name: annexin A13	
	Calculated MW: 316 aa, 35 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Annexin A13 concentration was determined to be 46.6 ng/mL in Jurkat cell extract based on a 1.3 mg/mL extract load.



Sandwich ELISA standard curve of MP02684-1, Human Annexin A13 Recombinant Matched Antibody Pair - PBS only. 85955-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag18581. 85955-2-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL