

For Research Use Only

# SLC22A15 Recombinant antibody, PBS Only

Catalog Number: 85960-1-PBS



## Basic Information

Catalog Number:

85960-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14547

GenBank Accession Number:

BC026358

GeneID (NCBI):

55356

UNIPROT ID:

Q8IZD6

Full Name:

solute carrier family 22, member 15

Calculated MW:

535aa, 59 kDa; 547aa, 61 kDa

Observed MW:

61-66 kDa

Purification Method:

Protein A purification

CloneNo.:

250347G3

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85960-1-RR (SLC22A15 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85960-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85960-1-RR against Human SLC22A15 were performed. The affinity constant is 5.00 nM.