

For Research Use Only

Mouse CtSD Recombinant antibody, PBS Only (Capture)

Catalog Number: 85964-3-PBS



Basic Information

Catalog Number:

85964-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2986

GenBank Accession Number:

NM_009983.3

GeneID (NCBI):

13033

UNIPROT ID:

P18242

Full Name:

cathepsin D

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

250329C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Cathepsin D concentration was determined to be 470.8 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02204-1, MOUSE Cathepsin D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85964-3-PBS. Detection antibody: 85964-2-PBS. Standard: Eg2986. Range: 0.781-100 ng/mL.