

For Research Use Only

ELA2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85965-1-PBS



Basic Information

Catalog Number:

85965-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3192

GenBank Accession Number:

NM_001972.4

GeneID (NCBI):

1991

UNIPROT ID:

P08246

Full Name:

elastase 2, neutrophil

Calculated MW:

29KD

Purification Method:

Protein A purification

CloneNo.:

250346C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Neutrophil Elastase concentration was determined to be 125.3 ng/mL in human serum.



Cytometric bead array standard curve of MP02206-1, Neutrophil Elastase Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85965-2-PBS. Detection antibody: 85965-1-PBS. Standard: Eg3192. Range: 0.625-80 ng/mL.