

For Research Use Only

Neutrophil Elastase Recombinant monoclonal antibody, PBS Only

Catalog Number: 85965-5-PBS



Basic Information

Catalog Number:

85965-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG39624

GenBank Accession Number:

NM_001972.4

GeneID (NCBI):

1991

UNIPROT ID:

P08246

Full Name:

elastase 2, neutrophil

Calculated MW:

29KD

Observed MW:

30 kDa

Purification Method:

Protein A purification

CloneNo.:

253414H2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

ELA2 (Neutrophil elastase) is also named as ELAL, NE, HLE, HNE, PMN-E, SERP1, GE, ELANE and belongs to the peptidase S1 family. It is a 33-kDa enzyme with several isoforms that differ in their extent of glycosylation (PMID:12223222). Human neutrophil elastase (HNE) is a serine protease with potent proteolytic activity (3), which is reported to cause mucin secretion (degranulation) (PMID:12169572). It may contribute to the directed migration of leukocytes into flamed postcapillary venules in addition to contributing to the process of tissue emigration (PMID:9683424). The full length protein has a signal peptide, propeptide and three glycosylation sites.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



U-937 cells were subjected to SDS PAGE followed by western blot with 85965-5-RR (ELA2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85965-5-PBS in a different storage buffer formulation.



HL-60 cells were subjected to SDS PAGE followed by western blot with 85965-5-RR (ELA2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85965-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85965-5-RR against Human Neutrophil Elastase were performed. The affinity constant is 2.51 nM.