

For Research Use Only

ADAMTS13 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85974-2-PBS



Basic Information

Catalog Number:

85974-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2556

GenBank Accession Number:

NM_139025

GeneID (NCBI):

11093

UNIPROT ID:

Q76LX8

Full Name:

ADAM metalloproteinase with
thrombospondin type 1 motif, 13

Calculated MW:

154 kDa

Purification Method:

Protein A purification

CloneNo.:

250498B6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ADAMTS13 concentration was determined to be 524.0 ng/mL in human serum.



Cytometric bead array standard curve of MP02235-1, ADAMTS13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85974-2-PBS. Detection antibody: 85974-1-PBS. Standard: Eg2556. Range: 0.781-100 ng/mL.