

For Research Use Only

PAPP A2 Recombinant antibody, PBS Only

Catalog Number: 85976-3-PBS



Basic Information

Catalog Number:

85976-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2812

GenBank Accession Number:

NM_020318.2

GeneID (NCBI):

60676

UNIPROT ID:

Q9BXP8-1

Full Name:

pappalysin 2

Calculated MW:

199 kDa

Observed MW:

245 kDa

Purification Method:

Protein A purification

CloneNo.:

250602C12

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 85976-3-RR (PAPPA2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85976-3-PBS in a different storage buffer formulation.