

For Research Use Only

AK1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85979-2-PBS



Basic Information

Catalog Number:

85979-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6890

GenBank Accession Number:

BC001116

GeneID (NCBI):

203

UNIPROT ID:

P00568

Full Name:

adenylate kinase 1

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

250557C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean AK1 concentration was determined to be 904.4 ng/mL in HeLa cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP02232-1, AK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85979-2-PBS. Detection antibody: 85979-1-PBS. Standard: Ag6890. Range: 0.781-100 ng/mL.