

For Research Use Only

ADAM10 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85989-4-PBS



Basic Information

Catalog Number:

85989-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23078

GenBank Accession Number:

BC126253

GeneID (NCBI):

102

UNIPROT ID:

O14672

Full Name:

ADAM metalloproteinase domain 10

Calculated MW:

748 aa, 84 kDa

Purification Method:

Protein A purification

CloneNo.:

250522A11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02389-1, Human ADAM10 Recombinant Matched Antibody Pair - PBS only. 85989-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag23078. 85989-4-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



Serum of six individual healthy human donors was measured. The ADAM10 concentration of detected samples was determined to be 698.4 pg/mL with a range of 232.3-2,831.4 pg/mL.