

For Research Use Only

RBM7 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85991-1-PBS



Basic Information

Catalog Number:

85991-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11106

GenBank Accession Number:

BC034381

GeneID (NCBI):

10179

UNIPROT ID:

Q9Y580

Full Name:

RNA binding motif protein 7

Calculated MW:

266 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

250586B1

Applications

Tested Applications:

IF/ICC, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02227-2, Human RBM7 Recombinant Matched Antibody Pair - PBS only. 85991-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag11106. 85991-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean RBM7 concentration was determined to be 3,857.0 pg/mL in HepG2 cell extract based on a 1.8 mg/mL extract load, 10,443.6 pg/mL in HeLa cell extract based on a 1.3 mg/mL extract load and 65,958.2 pg/mL in A549 cell extract based on a 3.1 mg/mL extract load.