

For Research Use Only

RBM7 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85991-2-PBS



Basic Information

Catalog Number:

85991-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11106

GenBank Accession Number:

BC034381

GeneID (NCBI):

10179

UNIPROT ID:

Q9Y580

Full Name:

RNA binding motif protein 7

Calculated MW:

266 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

250586E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean RBM7 concentration was determined to be 63.8 ng/mL in Hela cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP02227-1, RBM7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85991-3-PBS. Detection antibody: 85991-2-PBS. Standard: Ag11106. Range: 1.563-200 ng/mL