

For Research Use Only

S100A7 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85993-3-PBS



Basic Information

Catalog Number:

85993-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3738

GenBank Accession Number:

BC034687

GeneID (NCBI):

6278

UNIPROT ID:

P31151

Full Name:

S100 calcium binding protein A7

Calculated MW:

101 aa, 11 kDa

Purification Method:

Protein A purification

CloneNo.:

250425B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean S100A7 concentration was determined to be 242.9 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02225-1, S100A7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85993-3-PBS. Detection antibody: 85993-2-PBS. Standard: Ag3738. Range: 0.625-80 ng/mL.