

For Research Use Only

CFP Recombinant antibody, PBS Only (Detector)

Catalog Number: 85999-2-PBS



Basic Information

| | | |
|-------------------------------------|---|--|
| Catalog Number: 85999-2-PBS | GenBank Accession Number: BC015756 | Purification Method: Protein A purification |
| Concentration: 1 mg/ml | GeneID (NCBI): 5199 | CloneNo.: 250654E5 |
| Source: Rabbit | UNIPROT ID: P27918 | |
| Isotype: IgG | Full Name: complement factor properdin | |
| Immunogen Catalog Number: EG2219 | Calculated MW: 469 aa, 51 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Properdin/PFC concentration was determined to be 8.6 ug/mL in human serum.



Cytometric bead array standard curve of MP02228-1, Properdin/PFC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85999-3-PBS. Detection antibody: 85999-2-PBS. Standard: Eg2219. Range: 0.313-40 ng/mL.