

For Research Use Only

# CFP Recombinant antibody, PBS Only (Capture)

Catalog Number: 85999-3-PBS



## Basic Information

Catalog Number: 85999-3-PBS	GenBank Accession Number: BC015756	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 5199	CloneNo.: 250654E4
Source: Rabbit	UNIPROT ID: P27918	
Isotype: IgG	Full Name: complement factor properdin	
Immunogen Catalog Number: EG2219	Calculated MW: 469 aa, 51 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean Properdin/PFC concentration was determined to be 8.6 ug/mL in human serum.



Cytometric bead array standard curve of MP02228-1, Properdin/PFC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85999-3-PBS. Detection antibody: 85999-2-PBS. Standard: Eg2219. Range: 0.313-40 ng/mL.