

For Research Use Only

Mouse Ccl9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86002-1-PBS



Basic Information

Catalog Number:

86002-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3336

GenBank Accession Number:

NM_011338.2

GeneID (NCBI):

20308

UNIPROT ID:

Q3U9T8

Full Name:

chemokine (C-C motif) ligand 9

Calculated MW:

14KD

Purification Method:

Protein A purification

CloneNo.:

243103F7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CCL9/10/MIP-1 gamma concentration was determined to be 174.4 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02230-1, MOUSE CCL9/10/MIP-1 gamma Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86002-2-PBS. Detection antibody: 86002-1-PBS. Standard: Eg3336. Range: 1.563-100 ng/mL.