

For Research Use Only

Mouse Hpx Recombinant antibody, PBS Only (Capture)

Catalog Number: 86003-2-PBS



Basic Information

Catalog Number:

86003-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2829

GenBank Accession Number:

NM_017371.2

GeneID (NCBI):

15458

UNIPROT ID:

Q91X72

Full Name:

hemopexin

Calculated MW:

51 kDa

Purification Method:

Protein A purification

CloneNo.:

250600B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Hemopexin concentration was determined to be 110.9 ng/mL in Mouse brain tissue extract based on a 7.0 mg/mL extract load.



Cytometric bead array standard curve of MP02229-1, MOUSE Hemopexin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86003-2-PBS. Detection antibody: 86003-1-PBS. Standard: Eg2829. Range: 0.625-80 ng/mL.