

For Research Use Only

Mouse Hpx Recombinant antibody, PBS Only (Capture)

Catalog Number:86003-3-PBS



Basic Information

Catalog Number:

86003-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2829

GenBank Accession Number:

NM_017371.2

GeneID (NCBI):

15458

UNIPROT ID:

Q91X72

Full Name:

hemopexin

Calculated MW:

51 kDa

Purification Method:

Protein A purification

CloneNo.:

250600D5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02229-2, Mouse Hemopexin Recombinant Matched Antibody Pair - PBS only. 86003-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2829. 86003-1-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL



The mean Hemopexin concentration was determined to be 612.49 ng/mL in mouse brain tissue extract based on a 6.6 mg/mL extract load.