

For Research Use Only

HMGCR Recombinant antibody, PBS Only (Capture)

Catalog Number: 86005-2-PBS



Basic Information

Catalog Number:

86005-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4444

GenBank Accession Number:

BC033692

GeneID (NCBI):

3156

UNIPROT ID:

P04035

Full Name:

3-hydroxy-3-methylglutaryl-Coenzyme A reductase

Calculated MW:

888 aa, 95 kDa

Purification Method:

Protein A purification

CloneNo.:

243065G8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02239-1, Human HMGCRC Recombinant Matched Antibody Pair - PBS only. 86005-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag4444. 86005-1-PBS was HRP conjugated as the detection antibody. Range: 12.5-800 ng/mL



The mean HMGCRC concentration was determined to be 805.26 ng/mL in HeLa cell extract based on a 3.0 mg/mL extract load.