

For Research Use Only

ANGPTL4 Recombinant antibody

Catalog Number: 86006-4-RR



Basic Information

Catalog Number:

86006-4-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2770

GenBank Accession Number:

NM_020581.2

GeneID (NCBI):

57875

UNIPROT ID:

Q9Z1P8

Full Name:

angiopoietin-like 4

Calculated MW:

46 kDa

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

250335F5

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Positive Controls:

WB : rat liver tissue, mouse liver tissue, mouse brain tissue

Background Information

ANGPTL4 was first identified in adipose tissue and later found to be widely expressed in several tissues including the liver, skeletal muscle, small intestine, placenta, heart, and kidney. ANGPTL4 expression was first identified to be regulated by peroxisome proliferator-activator receptors (PPARs), PPAR α and PPAR γ .

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86006-4-RR (Angptl4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.