

For Research Use Only

# CCL27 Recombinant antibody

Catalog Number: 86007-2-RR



## Basic Information

Catalog Number:

86007-2-RR

Concentration:

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2461

GenBank Accession Number:

NM\_006664

GeneID (NCBI):

10850

UNIPROT ID:

Q9Y4X3

Full Name:

chemokine (C-C motif) ligand 27

Calculated MW:

13kd

Purification Method:

Protein A purification

CloneNo.:

242991B1

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

