

For Research Use Only

RBM47 Recombinant antibody, PBS Only

Catalog Number: 86009-1-PBS



Basic Information

Catalog Number: 86009-1-PBS	GenBank Accession Number: BC126261	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 54502	CloneNo.: 250436H8
Source: Rabbit	UNIPROT ID: A0AV96	
Isotype: IgG	Full Name: RNA binding motif protein 47	
Immunogen Catalog Number: AG20984	Calculated MW: 593 aa, 64 kDa	
	Observed MW: 64 kDa	

Applications

Tested Applications:
WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

RBM47, also known as RNA-binding motif protein 47 or RNA-binding protein 47, is a vertebrate-conserved RBP containing three RNA-recognition motifs (RRMs) that perform multifaceted roles in RNA editing, early embryonic development, and cancer. Recent studies suggest that RBM47 functions as a tumor suppressor. It has been shown to suppress the progression of breast and colon cancers. RBM47 represents a molecular switch of cell fate decisions that functions as a regulator of the p53/p21-signaling axis. (PMID: 34499322, PMID: 31511650)

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 86009-1-RR (RBM47 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86009-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 86009-1-RR (RBM47 antibody) at dilution of 1:500 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 86009-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 86009-1-RR (RBM47 antibody) at dilution of 1:500 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 86009-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A549 cells using RBM47 antibody (86009-1-RR, Clone: 250436H8) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 86009-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86009-1-RR against Human RBM47 were performed. The affinity constant is 1.02 nM.