

For Research Use Only

COPZ1 Recombinant antibody

Catalog Number: 86014-1-RR



Basic Information

Catalog Number:

86014-1-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14234

GenBank Accession Number:

BC002849

GeneID (NCBI):

22818

UNIPROT ID:

P61923

Full Name:

coatamer protein complex, subunit
zeta 1

Calculated MW:

177 aa, 20 kDa

Observed MW:

15-20 kDa

Purification Method:

Protein A purification

CloneNo.:

250388F3

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, monkey

Positive Controls:

WB : HeLa cells, RAW 264.7 cells, COS-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86014-1-RR (COPZ1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.