

For Research Use Only

# CD74 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86023-1-PBS



## Basic Information

Catalog Number:

86023-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1945

GenBank Accession Number:

NM\_001025159.3

GeneID (NCBI):

972

UNIPROT ID:

P04233-2

Full Name:

CD74 molecule, major histocompatibility complex, class II invariant chain

Calculated MW:

26kDa

Purification Method:

Protein A purification

CloneNo.:

250520A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean CD74 concentration was determined to be 28.6 ng/mL in human serum.



Cytometric bead array standard curve of MP02221-1, CD74 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86023-2-PBS. Detection antibody: 86023-1-PBS. Standard: Eg1945. Range: 1.563-200 ng/mL