

For Research Use Only

# MIA Recombinant antibody, PBS Only (Detector)

Catalog Number: 86040-1-PBS



## Basic Information

Catalog Number: 86040-1-PBS	GenBank Accession Number: NM_001202553.1	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 8190	CloneNo.: 250693B12
Source: Rabbit	UNIPROT ID: Q16674-1	
Isotype: IgG	Full Name: melanoma inhibitory activity	
Immunogen Catalog Number: EG2876	Calculated MW: 15 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean MIA concentration was determined to be 6.1 ng/mL in human serum.



Cytometric bead array standard curve of MP02259-1, MIA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86040-2-PBS. Detection antibody: 86040-1-PBS. Standard: Eg2876. Range: 0.781-100 ng/mL