

For Research Use Only

Mouse Lcat Recombinant antibody, PBS Only (Capture)

Catalog Number: 86041-2-PBS



Basic Information

Catalog Number: 86041-2-PBS	GenBank Accession Number: NM_008490.2	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 16816	CloneNo.: 250320A5
Source: Rabbit	UNIPROT ID: P16301	
Isotype: IgG	Full Name: lecithin cholesterol acyltransferase	
Immunogen Catalog Number: EG2756	Calculated MW: 50 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean LCAT concentration was determined to be 863.3 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02258-1, MOUSE LCAT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86041-2-PBS. Detection antibody: 86041-1-PBS. Standard: Eg2756. Range: 0.781-100 ng/mL