

For Research Use Only

# PSAP Recombinant antibody

Catalog Number: 86043-2-RR



## Basic Information

<b>Catalog Number:</b> 86043-2-RR	<b>GenBank Accession Number:</b> NM_011179.3	<b>Purification Method:</b> Protein A purification
<b>Concentration:</b> 1000 µg/ml	<b>GeneID (NCBI):</b> 19156	<b>CloneNo.:</b> 250312D3
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> NP_035309.3	<b>Recommended Dilutions:</b> WB: 1:2000-1:10000 IF/ICC: 1:200-1:800
<b>Isotype:</b> IgG	<b>Full Name:</b> prosaposin	
<b>Immunogen Catalog Number:</b> EG3155	<b>Calculated MW:</b> 61 kDa <b>Observed MW:</b> 65-70 kDa	

## Applications

<b>Tested Applications:</b> WB, IF/ICC, ELISA	<b>Positive Controls:</b> WB : NIH/3T3 cells, RAW 264.7 cells, mouse cerebellum tissue, mouse testis tissue IF/ICC : NIH/3T3 cells,
<b>Species Specificity:</b> mouse	

## Background Information

Prosaposin (PSAP) is also named as Sulfated glycoprotein 1 (SGP-1). PSAP is a secreted neurotrophic factor and a lysosomal protein serving as the precursor of saposins (A-D), which are cofactors for hydrolases of sphingolipids (PMID: 25130661, PMID: 31220974). PSAP modulates glycosphingolipid metabolism and variants have been linked to Parkinson's disease (PD) (PMID: 37726325). As a dual-function protein, PSAP is a lysosome-associated protein that participates in a variety of cellular processes. In the lysosome, PSAP is processed to activate enzymes that degrade lipids. PSAP proteins located extracellularly are involved in cancer progression, such as proliferation and tumor death suppression signaling (PMID: 39405903). Prosaposin (PSAP) has two forms: a precursor and a secreted form. The secreted form has neurotrophic, myelinotrophic, and myotrophic properties. The precursor form is a precursor protein of saposins A-D (PMID: 33242172). Both human and mouse prosaposin are 65-70 kDa proteins cleaved in the lysosomes into four 10-15 kDa proteins designated saposins A, B, C and D (PMID: 8573994). Sulfated glycoprotein-1 (prosaposin) exists in 2 forms: a 65 kDa form targeted to lysosomes and a post-translationally modified to the 70 kDa form secreted extracellularly (PMID: 8593668).

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3  
**Aliquoting is unnecessary for -20°C storage**

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86043-2-RR (Psap antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed NIH/3T3 cells using PSAP antibody (86043-2-RR, Clone: 250312D3 ) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).