

For Research Use Only

# Mouse Psap Recombinant antibody, PBS Only (Detector)

Catalog Number: 86043-4-PBS



## Basic Information

Catalog Number:

86043-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3155

GenBank Accession Number:

NM\_011179.3

GeneID (NCBI):

19156

UNIPROT ID:

NP\_035309.3

Full Name:

prosaposin

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

250312A8

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP02256-2, Mouse Psap Recombinant Matched Antibody Pair - PBS only. 86043-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3155. 86043-4-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL



Serum of sixteen mice was measured. The Psap concentration of detected samples was determined to be 396.8 ng/mL with a range of 2,614.3-9,786.9 ng/mL