

For Research Use Only

PPP1R13L Recombinant antibody, PBS Only (Capture)

Catalog Number: 86044-3-PBS



Basic Information

Catalog Number:

86044-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13273

GenBank Accession Number:

BC064913

GeneID (NCBI):

10848

UNIPROT ID:

Q8WUF5

Full Name:

protein phosphatase 1, regulatory
(inhibitor) subunit 13 like

Calculated MW:

89 kDa

Purification Method:

Protein A purification

CloneNo.:

250633C10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02334-1, Human iASPP Recombinant Matched Antibody Pair - PBS only. 86044-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag13273. 86044-1-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



The mean iASPP concentration was determined to be 32.51 ng/mL in PC-3 cell extract based on a 2.0 mg/mL extract load and 68.98 ng/mL in MCF-7 cell extract based on a 1.3 mg/mL extract load.