

For Research Use Only

IRF2BP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86048-3-PBS



Basic Information

Catalog Number:

86048-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4596

GenBank Accession Number:

BC038222

GeneID (NCBI):

26145

UNIPROT ID:

Q8IU81

Full Name:

interferon regulatory factor 2 binding protein 1

Calculated MW:

584 aa, 62 kDa

Purification Method:

Protein A purification

CloneNo.:

250630B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IRF2BP1 concentration was determined to be 27.1 ng/mL in Jurkat cell extract based on a 2.0 mg/mL extract load.



Cytometric bead array standard curve of MP02247-1, IRF2BP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86048-3-PBS. Detection antibody: 86048-2-PBS. Standard: Ag4596. Range: 0.781-100 ng/mL.