

For Research Use Only

Peripherin Recombinant antibody, PBS Only (Capture)

Catalog Number: 86057-2-PBS



Basic Information

Catalog Number:

86057-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11365

GenBank Accession Number:

BC032703

GeneID (NCBI):

5630

UNIPROT ID:

P41219

Full Name:

peripherin

Calculated MW:

470 aa, 54 kDa

Purification Method:

Protein A purification

CloneNo.:

250661F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Peripherin concentration was determined to be 14.5 ng/mL in COLO 320 cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP02242-1, Peripherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86057-2-PBS. Detection antibody: 86057-1-PBS. Standard: Ag11365. Range: 0.625-80 ng/mL.