

For Research Use Only

Peripherin Recombinant antibody, PBS Only (Capture)

Catalog Number: 86057-4-PBS



Basic Information

Catalog Number:

86057-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11365

GenBank Accession Number:

BC032703

GeneID (NCBI):

5630

UNIPROT ID:

P41219

Full Name:

peripherin

Calculated MW:

470 aa, 54 kDa

Purification Method:

Protein A purification

CloneNo.:

250661E8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02242-2, Human Peripherin Recombinant Matched Antibody Pair - PBS only. 86057-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag11365. 86057-1-PBS was HRP conjugated as the detection antibody. Range: 125-8000 pg/mL

The mean Peripherin concentration was determined to be 2.78 ng/mL in COLO 320 cell extract based on a 1.7 mg/mL extract load and 124.43 ng/mL in IMR-32 extract based on a 2.0 mg/mL extract load.