

For Research Use Only

# MAG/Siglec-4a Recombinant antibody

Catalog Number: 86073-1-RR



## Basic Information

Catalog Number:

86073-1-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3251

GenBank Accession Number:

NM\_017190.4

GeneID (NCBI):

29409

UNIPROT ID:

P07722-1

Full Name:

myelin-associated glycoprotein

Calculated MW:

69 kDa

Observed MW:

70-100 kDa

Purification Method:

Protein A purification

CloneNo.:

250740D11

Recommended Dilutions:

WB: 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Positive Controls:

WB : rat brain tissue, mouse brain tissue, rat spinal cord tissue, mouse spinal cord tissue, rat cerebellum tissue, mouse cerebellum tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86073-1-RR (MAG/Siglec-4a antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.