

For Research Use Only

PCNP Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86075-2-PBS



Basic Information

Catalog Number: 86075-2-PBS	GenBank Accession Number: BC013916	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 57092	CloneNo.: 250576E7
Isotype: IgG	UNIPROT ID: Q8WW12	
Immunogen Catalog Number: AG30951	Full Name: PEST proteolytic signal containing nuclear protein	
	Calculated MW: 19 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean PCNP concentration was determined to be 25.9 ng/mL in A2780 cell extract based on a 1.2 mg/mL extract load, 38.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02799-1, PCNP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86075-3-PBS. Detection antibody: 86075-2-PBS. Standard: Ag30951. Range: 0.938-60 ng/mL.



Bi-layer interferometry (BLI) kinetic assay of 86075-2-PBS against Human PCNP was performed. The affinity constant is 30.9 nM.