

For Research Use Only

SLC7A13 Recombinant antibody, PBS Only

Catalog Number: 86078-1-PBS



Basic Information

Catalog Number:

86078-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC125165

GeneID (NCBI):

157724

UNIPROT ID:

Q8TCU3

Full Name:

solute carrier family 7, (cationic amino acid transporter, y+ system) member 13

Calculated MW:

470 aa, 52 kDa

Observed MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

250062A4

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse

Background Information

SLC7A13 (also named AGT1 or XAT2) may enable L-amino acid transmembrane transporter activity. It may be involved in L-cystine transport; L-glutamate transmembrane transport; and aspartate transmembrane transport.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86078-1-RR (SLC7A13 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86078-1-PBS in a different storage buffer formulation.



IP result of anti-SLC7A13 (IP:86078-1-RR, 4ug; Detection:86078-1-RR 1:1000) with mouse kidney tissue lysate 3400 ug. This data was developed using the same antibody clone with 86078-1-PBS in a different storage buffer formulation.