

For Research Use Only

ASGR2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86079-2-PBS



Basic Information

Catalog Number:

86079-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3052

GenBank Accession Number:

N/A

GeneID (NCBI):

433

UNIPROT ID:

P07307-1

Full Name:

asialoglycoprotein receptor 2

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

250689A4

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02330-1, Human ASGR2 Recombinant Matched Antibody Pair - PBS only. 86079-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3052. 86079-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.



Serum of sixteen individual healthy human donors was measured. The ASGR2 concentration of detected samples was determined to be 4.46 ng/mL with a range of ND-18.65 ng/mL.