

For Research Use Only

ASGR2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86079-3-PBS



Basic Information

Catalog Number: 86079-3-PBS	GenBank Accession Number: N/A	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 433	CloneNo.: 250689H4
Source: Rabbit	UNIPROT ID: P07307-1	
Isotype: IgG	Full Name: asialoglycoprotein receptor 2	
Immunogen Catalog Number: EG3052	Calculated MW: 35 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02330-1, Human ASGR2 Recombinant Matched Antibody Pair - PBS only. 86079-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3052. 86079-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.



Serum of sixteen individual healthy human donors was measured. The ASGR2 concentration of detected samples was determined to be 4.46 ng/mL with a range of ND-18.65 ng/mL.