

For Research Use Only

Mouse Ccl8 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86081-2-PBS



Basic Information

Catalog Number:

86081-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2771

GenBank Accession Number:

NM_021443.3

GeneID (NCBI):

20307

UNIPROT ID:

Q9Z121

Full Name:

chemokine (C-C motif) ligand 8

Calculated MW:

11 kDa

Purification Method:

Protein A purification

CloneNo.:

250174H8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Ccl8 concentration was determined to be 1317.1 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02245-1, MOUSE Ccl8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86081-2-PBS. Detection antibody: 86081-1-PBS. Standard: Eg2771. Range: 0.156-20 ng/mL.