

For Research Use Only

MGLL Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86087-3-PBS



Basic Information

Catalog Number:

86087-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6970

GenBank Accession Number:

BC000551

GeneID (NCBI):

11343

UNIPROT ID:

A0A0C4DFN3

Full Name:

monoglyceride lipase

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

250509F8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MGLL concentration was determined to be 17.74 ng/mL in L02 cell extract based on a 3.90 mg/mL extract load and 5.47 ng/mL in SMMC-7721 cell extract based on a 2.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02391-1, Human MGLL Recombinant Matched Antibody Pair - PBS only. 86087-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag6970. 86087-3-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL