

For Research Use Only

PDS5B Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:86088-1-PBS



Basic Information

Catalog Number: 86088-1-PBS	GenBank Accession Number: NM_015032	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 23047	CloneNo.: 250461C9
Isotype: IgG	UNIPROT ID: Q9NTI5	
Immunogen Catalog Number: AG27977	Full Name: PDS5, regulator of cohesion maintenance, homolog B (S. cerevisiae)	
	Calculated MW: 165 aa	

Applications

Tested Applications:
IF/ICC, Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PDS5B concentration was determined to be 61.70 ng/mL in SKOV-3 cell extract based on a 1.20 mg/mL extract load and 55.20 ng/mL in HeLa cell extract based on a 0.90 mg/mL extract load.



Sandwich ELISA standard curve of MP02392-1, Human PDS5B Recombinant Matched Antibody Pair - PBS only. 86088-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag27977. 86088-3-PBS was HRP conjugated as the detection antibody. Range: 0.78-50 ng/mL