

For Research Use Only



Mouse FAM3D Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86090-4-PBS

Basic Information

Catalog Number: 86090-4-PBS	GenBank Accession Number: NM_146050.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 18300	CloneNo.: 250561F6
Isotype: IgG	UNIPROT ID: P97805	
Immunogen Catalog Number: EG2694	Full Name: oncoprotein induced transcript 1	
	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FAM3D concentration was determined to be 2291.8 ng/mL in Mouse large intestine tissue extract based on a 4.5 mg/mL extract load.



Cytometric bead array standard curve of MP02761-1, MOUSE FAM3D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86090-4-PBS. Detection antibody: 86090-1-PBS. Standard: Eg2694. Range: 0.312-40 ng/mL.