

For Research Use Only

PARP10 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86096-2-PBS



Basic Information

Catalog Number:

86096-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23457

GenBank Accession Number:

BC136542

GeneID (NCBI):

84875

UNIPROT ID:

Q53GL7

Full Name:

poly (ADP-ribose) polymerase family, member 10

Purification Method:

Protein A purification

CloneNo.:

25045483

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PARP10 concentration was determined to be 340.3 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 106.3 ng/mL in HL-60 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02820-1, PARP10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86096-3-PBS. Detection antibody: 86096-2-PBS. Standard: Ag23457. Range: 12.5-200 ng/mL.