

For Research Use Only

# RARG Recombinant antibody, PBS Only (Capture)

Catalog Number: 86098-2-PBS



## Basic Information

Catalog Number: 86098-2-PBS	GenBank Accession Number: BC019098	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 5916	CloneNo.: 250399B2
Source: Rabbit	UNIPROT ID: P13631	
Isotype: IgG	Full Name: retinoic acid receptor, gamma	
Immunogen Catalog Number: AG1991	Calculated MW: 454 aa, 50 kDa	

## Applications

Tested Applications:  
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP02317-1, Human RARG Recombinant Matched Antibody Pair - PBS only. 86098-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag1991. 86098-1-PBS was HRP conjugated as the detection antibody. Range: 46.9-3000 pg/mL.



The mean RARG concentration was determined to be 3.4 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.