

For Research Use Only

SNRPB Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86099-2-PBS



Basic Information

Catalog Number: 86099-2-PBS	GenBank Accession Number: BC003530	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 6628	CloneNo.: 250577G12
Isotype: IgG	UNIPROT ID: P14678	
Immunogen Catalog Number: AG10298	Full Name: small nuclear ribonucleoprotein polypeptides B and B1	
	Calculated MW: 234aa,25 kDa; 285aa,30 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean SNRPB concentration was determined to be 760.2 ng/mL in K-562 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02787-1, SNRPB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86099-3-PBS. Detection antibody: 86099-2-PBS. Standard: Ag10298. Range: 6.25-200 ng/mL.