

For Research Use Only

# A1BG Recombinant antibody, PBS Only (Capture)

Catalog Number: 86112-2-PBS



## Basic Information

Catalog Number:

86112-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1511

GenBank Accession Number:

NM\_130786.4

GeneID (NCBI):

1

UNIPROT ID:

P04217-1

Full Name:

alpha-1-B glycoprotein

Calculated MW:

54kDa

Purification Method:

Protein A purification

CloneNo.:

250024E6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean A1BG concentration was determined to be 91.4 ug/mL in human serum.



Cytometric bead array standard curve of MP02261-1, A1BG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86112-2-PBS. Detection antibody: 86112-1-PBS. Standard: Eg1511. Range: 1.563-100 ng/mL.