

For Research Use Only

Mouse VISTA Recombinant antibody, PBS Only (Detector)

Catalog Number: 86114-1-PBS



Basic Information

| | | |
|-------------------------------------|---|--|
| Catalog Number: 86114-1-PBS | GenBank Accession Number: NM_001159572.1 | Purification Method: Protein A purification |
| Concentration: 1 mg/ml | GeneID (NCBI): 74048 | CloneNo.: 250122F8 |
| Source: Rabbit | UNIPROT ID: Q9D659 | |
| Isotype: IgG | Full Name: RIKEN cDNA 4632428N05 gene | |
| Immunogen Catalog Number: EG2809 | Calculated MW: 34 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean VISTA concentration was determined to be 8.8 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02260-1, MOUSE VISTA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86114-2-PBS. Detection antibody: 86114-1-PBS. Standard: Eg2809. Range: 1.563-100 ng/mL