

For Research Use Only

LPP Recombinant antibody, PBS Only (Capture)

Catalog Number: 86117-2-PBS



Basic Information

Catalog Number:

86117-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19205

GenBank Accession Number:

BC130584

GeneID (NCBI):

4026

UNIPROT ID:

Q93052

Full Name:

LIM domain containing preferred translocation partner in lipoma

Calculated MW:

612 aa, 66 kDa

Purification Method:

Protein A purification

CloneNo.:

250773F6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean LPP concentration was determined to be 30.5 ng/mL in HeLa cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02282-1, LPP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86117-2-PBS. Detection antibody: 86117-1-PBS. Standard: Ag19205. Range: 0.156-20 ng/mL.