

For Research Use Only

Rat CD8a Recombinant antibody, PBS Only (Detector)

Catalog Number: 86119-2-PBS



Basic Information

Catalog Number:

86119-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3839

GenBank Accession Number:

NM_031538.2

GeneID (NCBI):

24930

UNIPROT ID:

P07725

Full Name:

CD8a molecule

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

250721H5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD8a concentration was determined to be 74.8 ng/mL in Rat spleen tissue extract based on a 8.8 mg/mL extract load.



Cytometric bead array standard curve of MP02416-1, RAT CD8a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86119-3-PBS. Detection antibody: 86119-2-PBS. Standard: Eg3839. Range: 0.469-60 ng/mL.