

For Research Use Only

# PRPF6 Recombinant antibody, PBS Only

Catalog Number: 86123-1-PBS



## Basic Information

Catalog Number:

86123-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13817

GenBank Accession Number:

BC001666

GeneID (NCBI):

24148

UNIPROT ID:

O94906

Full Name:

PRP6 pre-mRNA processing factor 6  
homolog (S. cerevisiae)

Calculated MW:

941 aa, 107 kDa

Observed MW:

100-108 kDa

Purification Method:

Protein A purification

CloneNo.:

250706F2

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

PRPF6 is involved in pre-mRNA splicing as component of the U4/U6-U5 tri-snRNP complex, it can enhance dihydrotestosterone-induced transactivation activity of AR, as well as dexamethasone-induced transactivation activity of NR3C1, but does not affect estrogen-induced transactivation.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86123-1-RR (PRPF6 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86123-1-PBS in a different storage buffer formulation.