

For Research Use Only

Prothrombin/F2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86124-2-PBS



Basic Information

Catalog Number:

86124-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1337

GenBank Accession Number:

NM_000506.5

GeneID (NCBI):

2147

UNIPROT ID:

P00734

Full Name:

coagulation factor II (prothrombin)

Calculated MW:

70kDa

Purification Method:

Protein A purification

CloneNo.:

243089B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Prothrombin/F2 concentration was determined to be 986.4 ug/mL in human serum.



Cytometric bead array standard curve of MP02277-1, Prothrombin/F2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86124-2-PBS. Detection antibody: 86124-1-PBS. Standard: Eg1337. Range: 0.781-100 ng/mL