

For Research Use Only

Prothrombin/F2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86124-3-PBS



Basic Information

Catalog Number:

86124-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1337

GenBank Accession Number:

NM_000506.5

GeneID (NCBI):

2147

UNIPROT ID:

P00734

Full Name:

coagulation factor II (prothrombin)

Calculated MW:

70kDa

Purification Method:

Protein A purification

CloneNo.:

243089B10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02277-2, Human Coagulation Factor II Recombinant Matched Antibody Pair - PBS only. 86124-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1337. 86124-2-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL



Plasma of sixteen individual healthy human donors was measured. The Coagulation Factor II concentration of detected samples was determined to be 64.2 µg/mL with a range of 40.1-77.6 µg/mL.